TORAY

TORAYSULFONE™ TS Series

Hollow Fiber Dialyzer
Dialyseur à fibres creuses
Dializzatore a fibre cave
Kapillardialysator
Dializador de Fibra Hueca
Φίλτρο αιμοκάθαρσης κοίλων ινών
Holle vezel dialysator
Dialisador de Fibra Oca

Instructions for Use of TORAYSULFONE
Instructions d'utilisation du TORAYSULFONE
Istruzioni per l'impiego del TORAYSULFONE
Gebrauchsanweisung für TORAYSULFONE Dialysatoren
Instrucciones para el Uso de TORAYSULFONE
Οδηγίες χρήσης του TORAYSULFONE
Gebruiksaanwijzing voor TORAYSULFONE
Instruções de Uso do TORAYSULFONE

EC REP

Toray International Italy S.r.l.

Via Mecenate 86, 20138 Milan, ITALY

Exporter:

Toray Medical Co., Ltd.

4-1, Nihonbashi-Honcho 2-chome, Chuo-ku, Tokyo 103-0023, JAPAN

Manufacturer:

Toray Industries, Inc.

1-1, Nihonbashi-Muromachi 2-chome, Chuo-ku, Tokyo 103-8666, JAPAN

ANVISA Registro:80975700003

Importador:

GSI CREOS BRASIL LTDA

Av. Paulista, 2001 room 214-216, São Paulo city 01311-931 Brazil

Exportador:

Toray Medical Co., Ltd.

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Reference code of issue date: {2403} ML

Instructions for Use of TORAYSULFONE TS Series

Read these instructions carefully before using the "TORAYSULFONE" TS series dialyzer.

I. INTENDED PURPOSE

"TORAYSULFONE" dialyzers are artificial kidney (hollow fiber dialyzer) which removes uremic toxins and water from the blood through hemodialysis to ameliorate conditions of a patient with renal insufficiency.

II. INDICATIONS

"TORAYSULFONE" dialyzers are indicated for SINGLE USE in acute or chronic hemodialysis. These dialyzers should be used only on the direction of a physician and used by trained medical personnel or by trained patient/caregiver authorized by a physician.

III. CONTRAINDICATIONS

Special contraindications for "TORAYSULFONE" are unknown. Generally, the contraindications for hemodialysis are applicable.

IV. STERILIZATION

"TORAYSULFONE" dialyzers are gamma-ray sterilized and nonpyrogenic (Blood side).

V. PERFORMANCE AND SPECIFICATIONS

Refer to the attached technical data.

VI. ADVERSE EVENTS

- 1) Patients should be carefully monitored during and after dialysis. In particular, patients with a history of allergy and hypersensitivity, patients who have experienced blood pressure decreasing by hemodialysis, patients whose immune function is accelerating, and patients using "TORAYSULFONE" for the first time should be carefully monitored. Appropriate treatments, including discontinuation of dialysis, should be taken, if the symptoms or signs of adverse reactions are observed.
- 2) Disequilibrium syndrome may appear especially during the introductory period of dialysis treatment, when applying to patients with low body weight, elderly patients, or patients in need of rapid water removal in a short period of time, or when changing dialyzer to one which is more efficient or has a larger effective surface area. Appropriate treatments, including discontinuation of dialysis, should be taken if the following symptoms or signs of disequilibrium syndrome are observed: headache, nausea, vomiting, consciousness disturbed, blood pressure decreased or convulsion.
- 3) Loss of proteins such as albumin may increase depending on the treatment conditions. Observe the patient's condition and serum albumin values.
- 4) Appropriate treatments, including discontinuation of dialysis, should be taken, if the following symptoms, which may occur occasionally in usual dialysis, are observed: dyspnea, chest pain, blood pressure decreased, blood pressure increased, loss of consciousness, hypoxemia, shock, palpitation, tachycardia, allergy, hypersensitivity, anaphylactoid reactions, white blood cell decreased, platelet decreased, eosinophilia, pruritus, edema (eyelid, throat, intraoral, face, etc.), malaise, fatigability, headache, vertigo, back pain, abdominal pain, lumbar pain, diarrhea, queasy, vomiting, chest discomfort, discomfort, eczema, rash, redness, urticaria, flushed face, pyrexia, rigors, abnormal sweating, muscle cramp, cough, hoarseness, tinnitus, yawning, abnormal vision, taste peculiar, strange smell sensation or hemolysis.
- 5) Improper handling or storage during transportation or operation (mechanical shock, excessively high or freezing temperature, etc.) may increase the risk of blood leakage. When blood leakage is found, hemodialysis should be discontinued and appropriate actions should be taken. The attending physician should decide whether or not to return the blood in the extracorporeal circuit to the patient.
- 6) Other complications may occur caused by malfunction, a defect of dialysis machine or a procedural error in its use, such as blood loss, blood overheating, hemolysis, excessive ultrafiltration and electrolyte imbalance. Refer to the instructions for the machine to avoid such complications.
- 7) In case of a serious incident related to this dialyzer, it should be reported to the manufacturer and the competent authority in which the user and/or patient is established.

VII. WARNINGS AND PRECAUTIONS

1. General warnings and precautions

Consult the physician if any abnormalities are found during and/or after the use of the dialyzer.

1) Reuse

"TORAYSULFONE" is a medical device intended for SINGLE USE ONLY. Its reuse might cause product damages (blood leakage, performance changes, etc.) and also induce

patient injury (infection caused by microbiological contamination, adverse reactions caused by residual reprocessing agents and/or blood components, etc.).

2) Special care

Special care is recommended for elderly, pregnant or pediatric patients.

3) Rinse

"TORAYSULFONE" dialyzers are filled with sterile water. Before starting dialysis, the blood side of the dialyzers must be flushed with sterile physiological saline or dialysate, and the dialysate side must be flushed with dialysate to avoid hemolysis.

4) Dialysis machine

Since "TORAYSULFONE" dialyzers have a high ultrafiltration rate, they should be used in conjunction with dialysis machines equipped with an ultrafiltration rate controller.

5) Storage

"TORAYSULFONE" should be stored and transported at temperatures between 5°C (41°F) and 30°C (86°F). Exposure to severe environmental conditions (direct sunlight, excessive humidity, wetness, dryness, freezing temperature) should be avoided.

6) Handling

"TORAYSULFONE" is made from various plastic components. Excessive vibration or mechanical shock should be avoided during handling and operation. Avoid tapping the dialyzer header with a clamp or other metal object to remove air bubbles or to return the blood.

7) Chemicals

Strictly avoid contact between the device and alcohol (sterilant) or other organic solvents. Such contact may lead to distortions or cracks of the device.

8) Dialysate quality control

Follow the standards of each country. Use clean dialysate since dialysate may come into patient blood stream through the mechanisms of backfiltration and /or backdiffusion.

2. Precautions to be taken before dialysis

- 1) The dialyzer must not be used if it appears to be damaged or its package is torn.
- 2) The dialyzer must not be used if there appears to be a water leak such as, drops of water in the package, improper cap fitting over either blood or dialysate port or an accumulation of a large quantity of air in the blood or dialysate compartment.
- 3) Connect the Luer lock connectors of the blood line to the blood ports by inserting at a correct angle and twisting to ensure a tight fitting to avoid detachment of connection, saline or blood leakage, or breakage of the blood ports.
- 4) The dialyzer should be used just after opening the package, rinsing and priming.
- 5) Take special care not to allow air bubbles to remain in or enter the blood compartment during the rinsing and priming procedure.

3. Precautions during dialysis

- 1) Confirm that there is no fluid leakage at the blood connection and the dialysate connection sites.
- 2) Transmembrane pressure (TMP) must not exceed 66 kPa (500 mmHg). TMP higher than 66 kPa (500 mmHg) could cause blood leakage or disconnection of the blood lines from the dialyzer connectors.
- 3) The performance of the dialyzer could be diminished if the device is used below the recommended flow rate or in orientations other than indicated.
- 4) The administration method and dose of an anticoagulant such as heparin should be determined for each patient by a physician. Special cares, including monitoring of coagulation time during dialysis, should be taken for patients who have a tendency to hemorrhage or who have a coagulation disorder.
- 5) If air bubbles, a blood leak, clotting or hemolysis are observed during dialysis, appropriate treatments, including discontinuation of dialysis or replacing the dialyzer, should be taken promptly.
- 6) When the administration of drugs is intended before or during dialysis, consider various conditions such as route, timing and dose, because dialysis may remove them or modify their effects. Special cares are recommended for patients receiving angiotensin-converting-enzyme inhibitors.
- 7) When completing dialysis, air rinsing of remaining blood should not be employed. Return the patient's blood by rinsing with saline, returning as much of the patient's blood as possible by holding the dialyzer vertically and rotating it around its axis, if necessary.

VIII. DIALYSIS PROCEDURE

Use protective gloves and an aseptic technique to prevent contamination of patient's blood when connecting the blood line to the dialyzer, collecting blood samples, or returning the blood from the

extracorporeal circuit. The extracorporeal circuit with the connectors conforming to ISO8638 or ISO8637-2 is recommended for this dialyzer.

The following procedure is an example.

1. Preparation for Dialysis

- 1) Take the dialyzer out of the package and place it on the dialyzer holder.
- 2) Hang a container of sterile physiological saline (1000 mL or more) on an IV pole, and connect the arterial infusion line to the container.
- 3) Fill the arterial blood line with physiological saline, first the vascular access connector side and then the dialyzer connector side. Clamp both ends of the arterial line. (Fig. 1)
- 4) Connect the venous line to the venous port of the dialyzer. (Fig. 2)
- 5) Carefully connect the arterial line filled with physiological saline to the arterial port taking care so as not to allow air to enter the dialyzer. If some air bubbles are observed in the arterial port, remove them by squeezing the venous drip chamber. (Fig. 3)
- 6) Rinse the blood side of the dialyzer and the blood line with more than 1000 mL of physiological saline at about 100 mL/min. (Fig. 4)
- 7) Connect the dialysate lines to the dialysate ports so that the dialysate flows countercurrently to the blood flow. Remove any air bubbles and rinse the dialysate side with dialysate flowing at about 500 mL/min. (Fig. 5) Let the dialysate flow for at least 5 minutes. Avoid any ultrafiltration during this procedure.
- 8) Prime the dialyzer and blood lines with more than 500 mL of heparinized physiological saline, and confirm that there are no air bubbles remaining in the dialyzer by flushing repeatedly. If massive air bubbles appear from fibers, replace the dialyzer. (Fig. 6)
- 9) Confirm that the venous bubble trap is about 3/4 full.
- 10) Clamp the venous line near the distal end.

2. Dialysis

- 1) Flow dialysate after setting filtration rate at 0 mL/h if dialysate hasn't flowed.
- 2) Connect the arterial line to the patient's arterial cannula.
- 3) Open the clamps of both the arterial and venous lines and start the blood pump at 50-100 mL/min. When blood reaches the end of the venous line, stop the blood pump and connect the venous line to the venous cannula.
- 4) Gradually increase the blood flow rate and the ultrafiltration rate to the levels required for each patient.
- 5) As for the confirmation of dialysis machine's alarm function prior to dialysis treatment, follow the machine manufacturer's instructions.

3. Completion of Dialysis

- 1) Reduce pressures in both the blood and dialysate sides as much as possible and stop the blood pump after 5 minutes of perfusion without ultrafiltration.
- 2) Connect a container containing 300 mL or more of physiological saline to the infusion line.
- 3) Using gravity, let physiological saline flow into the arterial blood access side to return in-line blood. Clamp the end of the line and aseptically withdraw the arterial cannula.
- 4) Remove the arterial blood line from the blood pump.
- 5) Using gravity and physiological saline, start blood return from the dialyzer and the blood lines.
- 6) Clamp the end of the line and aseptically withdraw the venous cannula.
- 7) Discard the dialyzer, the lines and other disposables in the manner approved by your institution.

IX. WARRANTY AND LIMITATION OF LIABILITY

- 1) "TORAYSULFONE" is manufactured in accordance with its specifications and in compliance with the GMP regulatory guidelines. Toray will replace a dialyzer with manufacturing defects free of charge. When necessary, Toray will request return of the dialyzer itself and all related packaging in order to investigate the cause and to take corrective action, if applicable.
- 2) Toray is not responsible for any damages, irrespective of product damages, patient injury or any other problems caused by misuse, improper handling, operation and storage by the customer.
- 3) Toray is not responsible for any product damages, patient injury or any other problems caused by reuse of the dialyzer.

X. OTHER INFORMATION

The following information is available on request.

- 1) Details of the test methods
- 2) *In vivo* performance characteristics

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