TORAY ETRF

TE-12R

Instructions for Use of TORAY ETRF

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Instructions for Use of TORAY ETRF

I. INTRODUCTION

To ensure the safe and effective operation of this device, the information provided in this leaflet must be read carefully before use.

WARNING: TORAY ETRF is not a dialyzer and must not be used as a haemodialyzer. This device must only be used with dialysis machines manufactured by Toray Medical Co., Ltd. and Toray Medical (Qingdao) Co., Ltd.

II. GENERAL INFORMATION

Intended use: TORAY ETRF is a sterile device designed to be used with a haemodialysis system for filtration of dialysis fluid during haemodialysis to remove bacteria and endotoxins and produce purified dialysis fluid for the procedure.

Storage: Store at a temperature between 5°C (41°F) and 30°C (86°F). Exposure to severe environmental conditions, such as direct sunlight, excessive humidity, moisture, excessively dry conditions or freezing temperatures, should be avoided.

Sterilization: This device has been sterilized by gamma-ray irradiation.

III. WARNINGS AND PRECAUTIONS

1) General warnings and precautions

- This device should be used only under the directions of a physician.
- TORAY ETRF is made from various plastic components. Excessive vibration or mechanical shock should be avoided during handling and operation.
- Avoid freezing the device because the device is filled with reverse osmosis (RO) water.
- Feed-water quality must be the same as or better than the quality of the RO water.
- Remove the device from its outer package just before use. The device must not be used if the following defects are found:
 - The outer package is damaged or stained;
 - The device is damaged, or other abnormalities are observed;
 - Plugs are not in place, or leakage is noted;
 - The expiry date has passed.
- Connect the device so that dialysis fluid is filtered from within the hollow fiber to the outside (Figure 1). If the device is connected incorrectly, adequate water quality may not be obtained.
- Rinse the device sufficiently on installation into the dialysis machine.
- The device must be disinfected and rinsed before and after use to maintain its performance appropriately. Confirm that the device is free from disinfectants before reuse.
- If the device is damaged, or other abnormalities are observed during use, immediately stop using it and replace it with a new one regardless of the expiry date.
- Even if the device is within the expiry date, the performance of the device may decline due to the excess use. Therefore, replace devices immediately after the period of use which is stated in this instruction for use.
- When hot water is used as disinfectant for a device where the expiry date has been exceeded, the device is likely to be damaged.
- Use protective measures including a gown, mask and gloves when handling this device.
- Discard this device after use in the manner specified by your institution.

2) During use

- Install and replace the device carefully to preserve sterility, and take special care to prevent air bubbles from remaining in the device.
- Use the device only after confirming that there has been no damage to the hollow fibers using the leak-checking function of the dialysis machine.
- Periodically drain the remaining fluid from the filter inlet compartment of the device. If this drainage is not performed, water quality may be inadequate. We recommend that you drain at least 300 mL of the remaining fluid before use.


Figure 1 Installation Method

Connect the device so that dialysis fluid is filtered from within the hollow fiber to the outside. Periodically drain the remaining fluid from the filter inlet compartment of the device.

IV. USE OF THE DEVICE

1) Installation

Consult the manufacturer of the dialysis machine for information about installation of the device. To install the device, follow the instructions provided in the manual of the dialysis machine.

2) Disinfection

TORAY ETRF must only be disinfected using the disinfectants listed in Table 1. Disinfectants not listed below must never be used. Carefully follow the disinfecting instructions provided in the dialysis machine manual.

Table 1 Disinfectants for TORAY ETRF.

	Concentrations	Duration (minutes)
Sodium hypochlorite	0.05–0.1%	20–60
Acetic acid	1% or less	20–60
Peracetic acid	0.02–0.04%	20–60
Hot water* (85°C (185°F) or below)	—	20–60

* Citric acid can be used in the hot water. The citric acid concentration must not exceed 2%.

V. EXPIRY DATE AND PERIOD OF USE

1) Expiry date

The device must be used before the expiry date indicated on the product label and outer case.

2) Period of use

- Replace devices after 100 treatment sessions or three months after installation, whichever comes first, on condition that the chemical disinfectants listed in Table 1 have been used. However, if there has been concurrent use of hot water (as listed in Table 1), then replace devices after 80 treatment sessions or three months after installation, whichever comes first.
- The aforementioned period of use is the recommended period of use. As this period of use is influenced by the conditions of use, it is important that you examine the water quality within the filter outlet compartment of the device periodically to confirm its retentive performance, and that you replace devices immediately if contamination is suspected.

SPECIFICATIONS AND PERFORMANCE

Hollow Fiber	Material	Polysulfone
	Inside Diameter (μm)	200
	Membrane Thickness (μm)	60
	Effective Surface Area (m ²)	1.2
Spacer Yarn		Polyethylene terephthalate
Housing Material		Polycarbonate
Potting Material		Polyurethane
O-ring Material		Silicone rubber
Sterilization		γ-ray irradiation
Performance Criteria		
Endotoxin Retentive Performance*,\$		LRV# ≥ 3
Bacterial Retentive Performance*,\$\$		LRV# ≥ 8
Max. TMP (Transmembrane pressure) (kPa (mmHg))		66 (500)

* in accordance with the voluntary standard of the Japan Medical Devices Manufacturers Association

\$. defined as $\log_{10}[\text{concentration of endotoxin (ET) in challenge solution}/\text{concentration of ET in filtrate}]$

\$\$: defined as $\log_{10}[\text{number of organisms in challenge suspension}/\text{number of organisms in filtrate}]$

LRV: Logarithmic reduction value

Upon request, information about methods used to obtain performance data is available.

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Toray International Italy S.r.l.
Toray Industries, Inc.